
## **Infrared sterilizer S825/S825L**

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# I: Brief introduction

**This brand infrared sterilizer** adopts international advanced design concepts and manufacturing technology, with smooth and beautiful appearance, compact and firm structure, high versatility, simple operation and other advantages. It is a rare instrument of scientific research and experiment. This instrument is mainly used for high-temperature sterilization and disinfection of small items such as inoculation loops and inoculation needles. It can replace alcohol lamps and is convenient and quick.

## II: Products overview

## III: Products introduction

## IV: Technical Data

### Model

**S825**

**S825L**

|   |   |   |
|---|---|---|
| Name                                    | Infrared Sterilizer   |   |
| Picture                                 |  |  |
| Voltage                                 | AC220V 50HZ   |   |
| Power                                   | 150   | 240   |
| Max.temperature inside the heating hole | 825°C±50°C  |   |
| Heating temperature                     | Two tempertures available 400°C/850°C   |   |
| Max-diameter of sterilizer goods        | Φ14   | Φ35   |
| Length of heating zone                  | 150   | 100   |
| Warm-up time                            | 10-20   |   |
| Fuse                                    | 250V 3A 5*20  |   |
| Ambient temperature                     | 5-40°C  |   |
| Relative humidity                       | ≤90%  |   |
| Product Size                            | 175*130*185   | 147*130*200   |
| Net weight                              | 1.3   | 1.4   |

# IV: Operation guide

## /// Basic operating guide

### Application

- > It is mainly used for microbiological experiments in biological safety cabinets, purification workbenches, mobile vehicles and other environments.

### Applied area

- > Used for high temperature sterilization of small items such as inoculation loops and inoculation needles

### Using Steps


- > Place the instrument on a stable workbench, file the switch in the OFF position
- > Plug in
- > Turn the switch to the required position

### Instrument use environment

- > Indoor use
- > The working environment of the instrument is suitable for the temperature range of +5°C~+40°C
- > The applicable temperature range of the instrument is ≤80%
- > There is no vibration and airflow that affect performance around
- > There is no conductive dust, explosive gas or corrosive gas in the surrounding air
- > Keep a safety distance of 10cm before and after the instrument is working

### Safety Precautions

- > It is strictly prohibited to plug in the power connector and toggle the power button when the hands are wet with liquid!
- > It is strictly prohibited to unplug the power plug when the instrument is live!
- > It is strictly forbidden to maintain and clean the instrument when it is powered on!
- > It is strictly forbidden to install the instrument on an uneven, swaying and vibrating work surface!

[www.labbox.com](http://www.labbox.com)