



INFORMACIÓN GENERAL

Nombre del producto : Pinzas de plástico para esmerilados
Descripción : De POM (Polioximetileno). No rayan el vidrio

DATOS TÉCNICOS

referencia	tipo	unidades por ref.
CHCP-014-010	naranja, para bocas 14/23	10
CHCP-019-010	azul, para bocas 19/26	10
CHCP-024-010	verde, para bocas 24/29	10
CHCP-029-010	rojo, para bocas 29/32	10

EMBALAJE Y DATOS LOGÍSTICOS

referencia	vol (cm3)	kg	TARIC	GTIN
CHCP-014-010	0,442	0,03	39269097	08434868013759
CHCP-019-010	0,3	0,04	39269097	08434868013766
CHCP-024-010	0,384	0,05	39269097	08434868013773
CHCP-029-010	1,188	0,08	39269097	08434868013797

FOTO DEL PRODUCTO

GENERAL INFORMATION

Product name : Conical head clip, plastic

Description : Manufactured in POM (Polyoxymethylene). Will not scratch glass surfaces

TECHNICAL DATA

reference	type	pcs/pack
CHCP-014-010	orange, for 14/23 joints	10
CHCP-019-010	blue, for 19/26 joints	10
CHCP-024-010	green, for 24/29 joints	10
CHCP-029-010	red, for 29/32 joints	10

PACKAGING AND LOGISTICS

reference	vol (cm3)	kg	TARIC	GTIN
CHCP-014-010	0,442	0,03	39269097	08434868013759
CHCP-019-010	0,3	0,04	39269097	08434868013766
CHCP-024-010	0,384	0,05	39269097	08434868013773
CHCP-029-010	1,188	0,08	39269097	08434868013797

PRODUCT PHOTO

INFORMATIONS GÉNÉRALES

Nom produit : Pince de rodage en plastique**Description :** Fabriquée en POM (Polyoxyméthylène). Ne raye pas le verre

DONNÉES TECHNIQUES

référence	type	unités par ref.
CHCP-014-010	orange, pour col 14/23	10
CHCP-019-010	bleu, pour col 19/26	10
CHCP-024-010	vert, pour col 24/29	10
CHCP-029-010	rouge, pour col 29/32	10

EMBALLAGE ET LOGISTIQUE

référence	vol (cm3)	kg	TARIC	GTIN
CHCP-014-010	0,442	0,03	39269097	08434868013759
CHCP-019-010	0,3	0,04	39269097	08434868013766
CHCP-024-010	0,384	0,05	39269097	08434868013773
CHCP-029-010	1,188	0,08	39269097	08434868013797

PHOTO PRODUIT

INFORMAZIONE GENERALE

Nome del prodotto : Clips in plastica per giunti conici smerigliati**Descrizione :** In POM (poliossimetilene). Non raschiano il vetro

DATI TECNICI

referenza	tipo	unità per ref.
CHCP-014-010	arancione, per coni 14/23	10
CHCP-019-010	blu, per coni 19/26	10
CHCP-024-010	verde, per coni 24/29	10
CHCP-029-010	rosso, per coni 29/32	10

IMBALLAGGIO E DATI LOGISTICI

referenza	vol (cm3)	kg	TARIC	GTIN
CHCP-014-010	0,442	0,03	39269097	08434868013759
CHCP-019-010	0,3	0,04	39269097	08434868013766
CHCP-024-010	0,384	0,05	39269097	08434868013773
CHCP-029-010	1,188	0,08	39269097	08434868013797

FOTO DEL PRODOTTO

ALGEMENE INFORMATIE**Produktnaam :** Kunststof klem voor slijpstukken**Beschrijving :** Van POM (Polyoxymethyleen). Laat geen krassen achter op het glas**TECHNISCHE GEGEVENS**

referentie	type	stuks per ref.
CHCP-014-010	oranje, voor NS 14/23	10
CHCP-019-010	blauw, voor NS 19/26	10
CHCP-024-010	groen, voor NS 24/29	10
CHCP-029-010	rood, voor NS 29/32	10

VERPAKKING EN LOGISTIEKE GEGEVENS

Referentie	vol (cm3)	kg	TARIC	GTIN
CHCP-014-010	0,442	0,03	39269097	08434868013759
CHCP-019-010	0,3	0,04	39269097	08434868013766
CHCP-024-010	0,384	0,05	39269097	08434868013773
CHCP-029-010	1,188	0,08	39269097	08434868013797

PRODUKTFOTO



ALLGEMEINE INFORMATIONEN

Produktname : Schliffklemmen aus Plastik
Beschreibung : Aus POM (Polyoxymethylene). Zerkratzt Glas nicht.

TECHNISCHE DATEN

Artikelnummer	Typ	Stückzahl pro Artikel
CHCP-014-010	orange, für NS 14/23	10
CHCP-019-010	blau, für NS 19/26	10
CHCP-024-010	grün, für NS 24/29	10
CHCP-029-010	rot, für NS 29/32	10

VERPACKUNG UND LOGISTIKDATEN

Referenz	vol (cm3)	kg	TARIC	GTIN
CHCP-014-010	0,442	0,03	39269097	08434868013759
CHCP-019-010	0,3	0,04	39269097	08434868013766
CHCP-024-010	0,384	0,05	39269097	08434868013773
CHCP-029-010	1,188	0,08	39269097	08434868013797

PRODUKTFOTO