



INFORMACIÓN GENERAL

Nombre del producto : Junta de goma Guko para matraces Kitasatos

Descripción : De goma natural de color gris


DATOS TÉCNICOS

referencia	talla	d (mm)	d2 (mm)	a (mm)	b (mm)	unidades por ref.
GUKO-001-001	1	21	12	18	2.5	1
GUKO-002-001	2	27	17	20	3	1
GUKO-003-001	3	33	21	24	3	1
GUKO-004-001	4	41	27,5	27	4	1
GUKO-005-001	5	53	33	34	5	1
GUKO-006-001	6	68	48	35	5.5	1
GUKO-007-001	7	78	58	35	6	1
GUKO-008-001	8	89	66	40	6.5	1

EMBALAJE Y DATOS LOGÍSTICOS

referencia	vol (cm3)	kg	TARIC	GTIN
GUKO-001-001	0,224	0,00	40169300	08434868025578
GUKO-002-001	0,336	0,00	40169300	08434868025592
GUKO-003-001	0,336	0,01	40169300	08434868025615
GUKO-004-001	0,448	0,01	40169300	08434868025639
GUKO-005-001	0,448	0,02	40169300	08434868025653
GUKO-006-001	2	0,04	40169300	08434868025677
GUKO-007-001	2	0,05	40169300	08434868025691
GUKO-008-001	2,25	0,07	40169300	08434868025714

FOTO DEL PRODUCTO



MATERIAL

MATERIAL Goma de color rojo



GENERAL INFORMATION

Product name : GUKO rubber gaskets, for filtering flask

Description : Made of natural rubber, grey


TECHNICAL DATA

reference	size	d (mm)	d2 (mm)	a (mm)	b (mm)	pcs/pack
GUKO-001-001	1	21	12	18	2.5	1
GUKO-002-001	2	27	17	20	3	1
GUKO-003-001	3	33	21	24	3	1
GUKO-004-001	4	41	27,5	27	4	1
GUKO-005-001	5	53	33	34	5	1
GUKO-006-001	6	68	48	35	5.5	1
GUKO-007-001	7	78	58	35	6	1
GUKO-008-001	8	89	66	40	6.5	1

PACKAGING AND LOGISTICS

reference	vol (cm3)	kg	TARIC	GTIN
GUKO-001-001	0,224	0,00	40169300	08434868025578
GUKO-002-001	0,336	0,00	40169300	08434868025592
GUKO-003-001	0,336	0,01	40169300	08434868025615
GUKO-004-001	0,448	0,01	40169300	08434868025639
GUKO-005-001	0,448	0,02	40169300	08434868025653
GUKO-006-001	2	0,04	40169300	08434868025677
GUKO-007-001	2	0,05	40169300	08434868025691
GUKO-008-001	2,25	0,07	40169300	08434868025714

PRODUCT PHOTO



MATERIAL

MATERIAL Red rubber

INFORMATIONS GÉNÉRALES

Nom produit : Joint en caoutchouc GUKO pour fiole à vide**Description :** En caoutchouc naturel, couleur grise


DONNÉES TECHNIQUES

référence	taille	d (mm)	d2 (mm)	a (mm)	b (mm)	unités par ref.
GUKO-001-001	1	21	12	18	2.5	1
GUKO-002-001	2	27	17	20	3	1
GUKO-003-001	3	33	21	24	3	1
GUKO-004-001	4	41	27,5	27	4	1
GUKO-005-001	5	53	33	34	5	1
GUKO-006-001	6	68	48	35	5.5	1
GUKO-007-001	7	78	58	35	6	1
GUKO-008-001	8	89	66	40	6.5	1

EMBALLAGE ET LOGISTIQUE

référence	vol (cm3)	kg	TARIC	GTIN
GUKO-001-001	0,224	0,00	40169300	08434868025578
GUKO-002-001	0,336	0,00	40169300	08434868025592
GUKO-003-001	0,336	0,01	40169300	08434868025615
GUKO-004-001	0,448	0,01	40169300	08434868025639
GUKO-005-001	0,448	0,02	40169300	08434868025653
GUKO-006-001	2	0,04	40169300	08434868025677
GUKO-007-001	2	0,05	40169300	08434868025691
GUKO-008-001	2,25	0,07	40169300	08434868025714

PHOTO PRODUIT



MATÉRIEL

MATÉRIEL Caoutchouc rouge

INFORMAZIONE GENERALE

Nome del prodotto : Guarnizione in gomma Guko per beute da vuoto**Descrizione :** In gomma naturale di colore grigio


DATI TECNICI

referenza	dimensione	d (mm)	d2 (mm)	a (mm)	b (mm)	unità per ref.
GUKO-001-001	1	21	12	18	2.5	1
GUKO-002-001	2	27	17	20	3	1
GUKO-003-001	3	33	21	24	3	1
GUKO-004-001	4	41	27,5	27	4	1
GUKO-005-001	5	53	33	34	5	1
GUKO-006-001	6	68	48	35	5.5	1
GUKO-007-001	7	78	58	35	6	1
GUKO-008-001	8	89	66	40	6.5	1

IMBALLAGGIO E DATI LOGISTICI

referenza	vol (cm3)	kg	TARIC	GTIN
GUKO-001-001	0,224	0,00	40169300	08434868025578
GUKO-002-001	0,336	0,00	40169300	08434868025592
GUKO-003-001	0,336	0,01	40169300	08434868025615
GUKO-004-001	0,448	0,01	40169300	08434868025639
GUKO-005-001	0,448	0,02	40169300	08434868025653
GUKO-006-001	2	0,04	40169300	08434868025677
GUKO-007-001	2	0,05	40169300	08434868025691
GUKO-008-001	2,25	0,07	40169300	08434868025714

FOTO DEL PRODOTTO



MATERIALE

MATERIALE Gomma di colore rosso

ALGEMENE INFORMATIE

Produktnaam : GUKO rubberen afdichting voor Kitasato-kolven**Beschrijving :** Van grijs natuurlijk rubber


TECHNISCHE GEGEVENS

referentie	maat	d (mm)	d2 (mm)	a (mm)	b (mm)	stuks per ref.
GUKO-001-001	1	21	12	18	2.5	1
GUKO-002-001	2	27	17	20	3	1
GUKO-003-001	3	33	21	24	3	1
GUKO-004-001	4	41	27,5	27	4	1
GUKO-005-001	5	53	33	34	5	1
GUKO-006-001	6	68	48	35	5.5	1
GUKO-007-001	7	78	58	35	6	1
GUKO-008-001	8	89	66	40	6.5	1

VERPAKKING EN LOGISTIEKE GEGEVENS

Referentie	vol (cm3)	kg	TARIC	GTIN
GUKO-001-001	0,224	0,00	40169300	08434868025578
GUKO-002-001	0,336	0,00	40169300	08434868025592
GUKO-003-001	0,336	0,01	40169300	08434868025615
GUKO-004-001	0,448	0,01	40169300	08434868025639
GUKO-005-001	0,448	0,02	40169300	08434868025653
GUKO-006-001	2	0,04	40169300	08434868025677
GUKO-007-001	2	0,05	40169300	08434868025691
GUKO-008-001	2,25	0,07	40169300	08434868025714

PRODUKTFOTO



MATERIAAL

MATERIAAL Rood rubber



ALLGEMEINE INFORMATIONEN

Produktname : Guko-Gummidichtungen für Saugflaschen
Beschreibung : Aus grauem Naturgummi.


TECHNISCHE DATEN

Artikelnummer	Größe	d (mm)	d2 (mm)	a (mm)	b (mm)	Stückzahlpro Artikel
GUKO-001-001	1	21	12	18	2.5	1
GUKO-002-001	2	27	17	20	3	1
GUKO-003-001	3	33	21	24	3	1
GUKO-004-001	4	41	27,5	27	4	1
GUKO-005-001	5	53	33	34	5	1
GUKO-006-001	6	68	48	35	5.5	1
GUKO-007-001	7	78	58	35	6	1
GUKO-008-001	8	89	66	40	6.5	1

VERPACKUNG UND LOGISTIKDATEN

Referenz	vol (cm3)	kg	TARIC	GTIN
GUKO-001-001	0,224	0,00	40169300	08434868025578
GUKO-002-001	0,336	0,00	40169300	08434868025592
GUKO-003-001	0,336	0,01	40169300	08434868025615
GUKO-004-001	0,448	0,01	40169300	08434868025639
GUKO-005-001	0,448	0,02	40169300	08434868025653
GUKO-006-001	2	0,04	40169300	08434868025677
GUKO-007-001	2	0,05	40169300	08434868025691
GUKO-008-001	2,25	0,07	40169300	08434868025714

PRODUKTFOTO



MATERIAL

MATERIAL Roter Gummi